

## FACT SHEET CURE TAPE

### MEDICAL DEVICE STATUS AND CERTIFICATION

Cure Tape is registered in Holland as a Class I Medical Device: registration number is NL/CA01/04-07434.

### WHAT IS THE CURE TAPE PRODUCT LINE?

The Cure Tape Product line is used by physiotherapists, nurses, osteopaths, sports massage therapists, occupational therapists and medical doctors to complement existing techniques in physical rehabilitation and sports injuries.

The basic methods of application (known as Medical Taping Concept or Neuromuscular Taping Method) is taught in a 15 hour course.

Cure Tape has the following product characteristics:

- a. **There is no chemical or pharmacological product delivered by the tape.**
- b. **The adhesive is alcohol free with hypo-allergenic backing, in wave-form, reducing the risk of skin irritation and allowing the skin to breathe.**
- c. **The tape is not latex but based on cotton fibre; this provides elasticity similar to that of the human skin (130%-140% self-stretching).**
- d. **The adhesive characteristics of the surface are activated by body heat.**
- e. **The tape is air and water-permeable, while still being water resistant; allowing for the tape to be worn for several days (4-5 days) without loss in the quality of its effect.**

### OBJECTIVES IN THE APPLICATION OF CURE TAPE:

- a. **To raise the skin to enhance both blood and lymphatic circulation.**
- b. **To support muscles, ligaments and tendons during rehabilitation or recovery from sports injuries.**

### PRODUCT LINE DESCRIPTION

	NAME	DESCRIPTION
1	Cure Tape Natural	Neuromuscular Tape 5,0 cm * 5,0 m-water resistant
2	Cure Tape Blue	Neuromuscular Tape 5,0 cm * 5,0 m-water resistant
3	Cure Tape Pink	Neuromuscular Tape 5,0 cm * 5,0 m-water resistant
4	Cure Tape Black	Neuromuscular Tape 5,0 cm * 5,0 m-water resistant
5	Cure Tape Orange	Neuromuscular Tape 5,0 cm * 5,0 m-water resistant
6	Cure Tape Natural - 7,5 cm	Neuromuscular Tape 7,5 cm * 5,0 m-water resistant

### SHELF LIFE OF CURE TAPE

The Cure Tape product shelf life is 30 months and should be stored in a cool area to preserve adhesive qualities.

### TAPING TECHNIQUES

While the neurological taping concept was originally developed in Japan and Korea in the 1970's, Cure Tape represents the "hardware" behind the Medical Taping concept. The "software" behind the Concept continues to evolve. The evolution has led to the development of the following taping techniques:

- Muscle Tape
- Facia Tape
- Ligament/Tendon Tape
- Corrective Techniques

- Space Technique
- Lymph Tape

These various techniques are taught in the following 10 to 15 hour modular courses scheduled for:

- **2011 – Existing Courses:** Basic Course; Sports Course; Advanced Course; Lymphatic Course; Pediatric Course and Neurological Course

For more information on courses, please see [www.aneid.com](http://www.aneid.com) or contact Mycology Research Laboratories (trading as Aneid UK ) at the following contact numbers: [Tel:\(0\)1582-485-209](tel:01582-485-209) or email [info@mycologyresearch.com](mailto:info@mycologyresearch.com) .